

ANALYTICAL BALANCE

0.001g / 0.01g

TOPLAB INDIA
"Serving Science with Solution"



The visual appearance of the product is blended into the simple concept, stylish and generous. Sliding glass wind screen provides big space, insightful vision and operates easily. TL series is the dual range.

TECHNICAL SPECIFICATION:

| MODEL | TL-103 | TL-203 | TL-303 | TL-106 |
|-------------------|---|--------|-----------------------|--------------------|
| Capacity (g) | 110 | 210 | 310 | 100/600 dual range |
| Sensor | Load Cell | | | Load Cell |
| Resolution (g) | 0.001 | | | 0.001/0.01 |
| Min weighing (g) | 0.010 | | | 0.010/0.04 |
| Stable Time | ≤ 3 S | | | |
| Display | LCD (white back light with black font) | | | |
| Pan Size | Φ 90mm | | Φ 90mm/ Φ 110mm | |
| Operation temp | 5-25°C | | | |
| Repeat ability | ±0.002g | | ±2d | |
| Liner | ±0.003g | | ±3d | |
| Draft shield size | 177mmx173mmx205mm | | | |
| Cal. weight | 100g | 200g | | 200g |
| Option | RS232C & RS485, Printer, Dual display, USB... | | | |
| Packing Size | Inside:350×265×345mm | | Outside:740×550×390mm | |
| G.weight | 4.0/16kg | | | |

ANALYTICAL BALANCE

